

NEW



Laminectomy Punches with **15 mm opening jaw** and semi-detachable shaft

Punches with 15 mm opening jaw for spinal vertebral operations. This longer opening jaw allows easier mobility of the instrument in a narrow operation field, thus allowing access to larger bones. Thereby the bone tissue can be eliminated step by step in an efficient way.



Product characteristics:

- Punches with 15 mm opening jaw
- Jaw 130° upwards cutting
- With semi-detachable shaft
- Without ejector
- From 2 to 5 mm
- Length 200 mm
- With flat footplate
- With or without ceramic coating

Width of bite	Length of shaft	 upwards cutting 130° without ceramic coating	 upwards cutting 130° with ceramic coating
2 mm	200 mm	▶ 45-1483	45-1483-C
3 mm	200 mm	45-1484	45-1484-C
4 mm	200 mm	45-1485	45-1485-C
5 mm	200 mm	45-1486	45-1486-C

Price information on request!