


In addition to our successful line „Laminectomy Punches with semi-detachable shaft“ we now would like to present an ideal complement for the neurosurgery.

## Rongeurs for Intervertebral Discs with semi-detachable shaft

The easy and fast opening and closing of the shaft allow the best possible cleaning and therefore it reduces the risk of infection.



Rongeurs for Intervertebral Discs CUSHING, LOVE-GRUENWALD, SPURLING type (24 Models)



Length of shaft: 180 / 200 mm  
Size of jaw: 2 - 5 mm  
Types: straight/curved upwards/curved downwards

Rongeurs for Intervertebral Discs CASPAR type (12 Models)



Length of shaft: 180 / 200 mm  
Size of jaw: 2 - 4 mm  
Types: straight/curved upwards

### Fax-Answer

Please send us: ..... catalogue(s) „Rongeurs for Intervertebral Discs No. 09/125“ - Cat. No. 100-0136

-----  
date

-----  
company's name

-----  
contact person